

Reizen Hypoglycemia Watch (Item 188221)

GETTING STARTED:

1. Insert Battery. Using a small screwdriver, slide battery drawer on side of watch open. Insert one CR2032 battery with Positive (+) side facing down

2. Put It On. Put the Hypoglycemia Watch on your wrist or ankle using the elastic strap/band. (Note: The two (2) metal sensors on the watch back must be in contact with the skin.)

3. Activate It. To turn the monitor on, press the On/Off button on the front of the monitor watch once. You should hear three (3) quick beeps indicating the monitor is "On".

4. Check Status. Press On/Off button to verify monitor is on. You will hear two quick beeps to confirm it is on.



ALARMS/ALERTS:

- **Alarm/Alert.** When the alarm is "On" (in active mode), it will sound four (4) continuous beeps and vibrate to alert the user to either the presence of perspiration or a drop in skin temperature.

- **Turn Alarm Off (Alarm Sounding).** To silence the alarm/alert while it is sounding, press On/Off button twice quickly (within one second). You will hear one long beep to indicate alert has been turned off.


- **Test Blood Sugar.** If an alarm sounds, test your blood sugar level using a blood glucose monitor to determine if hypoglycemia is present.

TURNING OFF:

- **Turn Off (Alarm Not Sounding).** Press and hold the On/Off button down for ten (10) seconds until you hear an alarm, then press On/Off button twice quickly within one second. You will hear one long beep indicating monitor is "Off".

RE-SETTING:

- **Reset Monitor Watch.** Remove battery and hold On/Off button down for ten (10) seconds. Replace battery and resume normal operation.



MaxiAIDS

PRODUCTS FOR INDEPENDENT LIVING

Call: 1-800-522-6294 • Fax: 631-752-0689
TTY: 631-752-0738 • Web: www.MaxiAids.com